

How to use the Zetasizer Nano ZS ?

- Switch on the computer and log in on the local session .\dls_malvern with zetasizer@MHMC as password
- Login to: <https://imxinstruments1.epfl.ch/cae/cae.html> with your Gaspar credentials to switch on the instrument
- For disposable cuvettes : Fill in to approx. 1 cm from the bottom
For the DTS1070, fill it upside down until the bend to avoid bubbles and then to the fill max limit
- Open the Zetasizer Software
- Click on Measure then Manual
- Select the measurement type (Size or Zeta potential)
- Under Material, either choose an existing material or create your own but **do not modify an existing material**
- Same with Dispersant
- Leave the General options as default
- Temperature : set the measurement temperature (max 70°C) and equilibration time
- Cell : Select right cuvette (DTS1070 for Zeta potential)
- Measurement : Leave as default
- Press Ok then press the Start button
- During the measurement, check the Multi-view tab to have an idea of the quality of the measure (cf. Chapter 4 of the User Manual)
- Also, at the bottom right, check the Attenuator used (Attn:)
 - Ideal range from 6 to 9.
 - If > 9 sample too dilute and < 6 sample too concentrated
- Once you're done with your measurements, **empty the cuvette content in the liquid waste container**, and if using the DTS1070, wash it with MQ water
- Close the software and turn off the instrument on the CAE