


EPFL - STI - CMi

ALPHA STEP 500

The Alpha-Step 500 equipment is used to characterize a sample surface profile by scanning it with a diamond tip sensor. The obtained plot is a representation of the scanned area cross-view.

It can measure steps height down to 50Å through a scanning length of maximum 10mm. The holder of this particular equipment is suitable wafers up to 100mm.

Do **NOT** put the tip in a hole or at any position away from the wafer. To measure steps higher than 500µm, always do it in descent (from the top of the feature):




- 1) Switch on the computer screen
- 2) Press the **Z-θ** key and check that the holder is in low position
- 3) Put your sample on the holder
- 4) Activate the vacuum clamping (*vac* position) and press the **LOAD** key
- 5) Keep the ↓ key pressed until the tip contacts the sample

NEVER MOVE THE HOLDER WHEN THE TIP IS IN CONTACT WITH THE SAMPLE

- 6) Press once the ↑ key to slightly lift the tip from the surface
- 7) Move to the place you want to measure using the “*trackball*” (use **F1** key to change the moving speed)
- 8) Press **F4** or **ESC** keys to return to the main menu

- 9) Choose CATALOG in the RECIPE menu. Either choose the general recipe or make your own (1 recipe per user). See appendices for recipe writing
- 10) Return to “camera view” using the **Z-θ** key
- 11) Press the **START** key to initiate the measure
- 12) Level adjustment of the profile plot
 - Press the **LEVEL** key
 - Using the 2 displayed red cursors, indicate two points on the profile plot that should be at the same level. Use the ← and → **keys** to move the cursors and the **space bar** to switch from one to another
 - Press **ENTER**. The profile plot level is adjusted



- 13) Measuring cursors
 - Use the two cursors to get the height between two points of the plot
 - Average the value targeted by each of the two cursors by splitting them using the ↑ and ↓ keys
 - Use the **space bar** to switch from one cursor to the other or to activate both
- 14) Profile zoom
 - Press **ZOOM** (a “zoom box” appears)
 - Use the ← and → **keys** or the “**trackball**” to move that box to the area you want to zoom in

- Press **F1** to anchor the box and use the ← and → *keys* or the “*trackball*” to adjust the “zoom box” size
- Press **F2** or **ZOOM** for seeing (Press **F3** to get back to original plot)
- Press **F4** to suppress the “zoom box”

15) Data saving

- Press **F4** (save)
- “DATA BASE SAVE IDs” window opens
- Enter # symbol into the first six ID fields (quick shortcut. Press **F1:Default**)
- In the “part ID” field enter the name as follows : #**nom**
- In the “Save data ?” field choose **SUMM + RAW**
- Press **ENTER** to confirm

16) Data exporting on floppy disk

- Press **F4** or **ESC** to return to main menu
- Choose DISK PATH in CONFIG menu and press **ENTER** to see DISK PATH window
- Select “Export to/Import from: C:” and press **F4**
- Select CATALOG in DATA menu and press **ENTER**
- In “Type” field enter “Any Data” and press **F2**
- Select the file and press **F3(export)**
- Press **N** (ASCII format)
- Type a file name and press **ENTER**

17) End of operation

- Go back to camera view by pressing the **Z-θ** key and then press **UNLOAD**
- Release the vacuum clamping and unload your sample
- Go back to main menu by pressing **F4** or **ESC**
- Switch off the screen (do not switch off the computer, only the screen)