

## *Product Information Brochure*



### *Shadow Mask Aligner*

Idonus Sarl  
Rouges-Terres 61  
2068 Hauterive / NE  
Switzerland

Phone : +41 32 724 44 40  
Fax : +41 32 724 44 42  
Email : [info@idonus.com](mailto:info@idonus.com)  
Web : [www.idonus.com](http://www.idonus.com)

## Shadow Mask Aligner

### General Information

Thin film deposition through a shadow mask is a powerful process for PVD layer structuring. It allows to replace one photolithographic and one etching step by precise masking of the substrate during deposition. A basic problem of the process is the alignment and temporary fixation of the shadow mask and the wafer.

Idonus provides a system for accurate and reproducible alignment of a shadow mask and a wafer with an integrated vibration-resistant clamping chuck. An alignment tolerance of 6  $\mu\text{m}$  was demonstrated with the help of a microscope. The chuck can directly be inserted in the vacuum chamber for PVD. After deposition the chuck is opened and the wafer and the mask are separated without damaging the wafer. The shadow mask can be reused.

Deposition chucks can be supplied for any substrate size and can be adapted to your requirements. We also provide precise and cost efficient shadow masks made of various materials like metal, silicon, glass or ceramic.




### Benefits of the Shadow Mask Aligner

- Precise and reproducible PVD layer structuring
- Vibration-resistant chuck for PVD
- Easy, damage free wafer clamping mechanism
- Simple handling during alignment due to vacuum clamping
- Ideal for different wafer sizes
- Customized chuck
- No installation
- Low running costs

# Shadow Mask Aligner

## Schema

innovations for microsystem fabrication



Principle section view of the complete system






Image of a 100 mm system with a chuck ready for alignment

## Shadow Mask Aligner

### Details: Chuck



Principle section view of the chuck with clamped wafer and mask. With the chuck any wafer and mask thickness can be clamped. The chuck exists for 100 / 150 and 200mm standard wafer size. Special sizes can be supplied.




Bottom (left) and top part of the chuck with vacuum connector. For the bottom part we also fabricate a version with opening.

## Shadow Mask Aligner

*Details: Vacuum unit*

innovations for microsystem fabrication



Bottom unit slid out to load bottom chuck with wafer.  
Top chuck ready to insert mask (mask not loaded in picture).



Pneumatic board for individual regulation of wafer,  
mask and chuck clamping during alignment.

# Shadow Mask Aligner

## Technical Specifications

Product Code	SMA
<b>Wafer and mask sizes</b> Diameter Thickness	100, 150 or 200 mm depending on equipment configuration independent
<b>Stage characteristics</b> Range in X and Y Range in Z Rotation Tilt	Scale step 10 µm - sensitivity 2 µm Scale step 10 µm - sensitivity 2 µm 360°endless Factory calibrated and manual adjustment by technician
<b>Alignment precision</b> Linear alignment Rotation alignment	After sputtering Typically 6 µm ; Max 10 µm 0.001 to 0.009 Deg
<b>Chuck</b> Wafer sizes Chuck material Chuck thickness	100 mm or 150 mm or 200 mm Aluminium with or without anodizing 20 mm depending on wafer and mask thickness
<b>Shadow mask</b> Material	Customized shadow masks can be provided Silicon, quartz, glass, metal, ceramic
<b>Dimension (mm)</b> Footprint (WxD)	500 x 250 mm <sup>2</sup>
<b>Installation</b> Need of	Microscope for visual control during alignment Vacuum

**Note:** System can be customized, please sent us your special requirements!