

# Dataplate User Manual

Version of 2020-08-04.

## 1. Introduction

This user manual explains how to operate the Dataplate hotplate for baking steps of SU-8 films.

## 2. Login

- Login on “Z01 DataPlate” with CAE on zone 01 accounting computer.


Z01 DataPlate - Programable Hotplate

## 3. Operation

- The hotplate uses its own controller with an extremely unintuitive ☹ language.
- Make sure the controller is turned on. If not, turn it on.
- The Dataplate can operate in two modes:
  - Fixed temperature mode
  - Programmed temperature mode, preferred for SU-8 bake steps.


**WARNING: Maximum hotplate temperature is 220°C.**

### 1) “Fixed temperature” mode:





- The display shows the actual temperature.
- Users can change the temperature (button 1) and timers (buttons 4 or 5) settings by using the right combination of key presses.

- Temperature: For 120°C setpoint, press:




- Timer: For 15 minutes, press:




- Pressing **button “1”** will start the **heating**. The status LED will turn ON.
- Make sure to wait at least 5 minutes for temperature stabilization.
- Open the cover.
- Place your wafers on the aluminium foil of the hotplate (maximum 5 wafers in a batch).
- Close the cover.
- Pressing **buttons “4” or “5”** will start the **timers**. The status LED will turn ON.
- When the timer is finished, turn off the hotplate.

### 2) “Programmed temperature” mode:

- The “Programmed temperature” mode uses the same concept of key presses to edit the temperature profile program, but with much more complexity.
- Turn on the hotplate.
- Enter the program edition by pressing:




- The “edit” status LED will turn ON. Step “00.” appears. **At this step, pressing “CLR” will erase the complete program.**




- Buttons “7” and “8” will cycle through the steps of the program.




- At any steps other than “00.”, pressing “CLR” will erase the actual step.
- Each step of the program can set one of the four following parameters:
  - button “1” = temperature [°C]
  - button “4” = timer [hh:mm]
  - button “5” = timer [mm:ss]
  - button “8” = ramp speed [°/hour]
- Each step will have one of two different “functions” which will be toggled by pressing several times (before the parameter keys) on:



- The function N° 1 is “SET”, the display symbol is “|\_”: **Immediately set target value and continue to the next step.**
- The function N° 2 = “WAIT UNTIL”, the display symbol is “|\_|”: **Continue to the next step when the target is reached.**
- Button “9” = HEATER OFF function
- When all program steps are entered correctly, users can out of the edition mode with button “6”. The status LED will turn OFF.



- The program is started by pressing “CLR” (outside of edition mode)




#### 4. Program example: SU-8 GM1075 PEB

- Please find in the table below the recommended parameters from the SU-8 runcard for the GM1075 PEB step.

Z1 / Accuplate Z1 / Dataplate	<u>Segment 1:</u> Ramp time: X min Temp.: Z °C; Step time: Y min  <u>Segment 2:</u> Ramp time: 90 min Temp.: 30°C; Step time: 0 min  <u>Segment 3 -&gt; 5:</u> NO INFLUENCE	X is the time for the ramp up, Y is the bake time.  [Type: X, Y] 1040: 20 min, 15 min 1050: 20 min, 20 min 1060: 30 min, 30 min 1070: 40 min, 40 min 1075: 50 min, 60 min  Z is the bake temperature: 90C for Si 85C for float glass
--	--	--

- Enter edition mode and erase the program:




- STEP “01.” The starting temperature is 20°C and the first ramp is 50 minutes to go to 90°C (Si wafer) → ramp speed 84°/hour. The following buttons sets the ramp speed.



Make sure the function is “|\_”.

- STEP “02.” The ramp speed is fixed, we can start heating until we reach 90°C. The following buttons sets the heating ramp until 90°C.



Make sure the function is “|\_|”.

- STEP “03.” We reached 90°C and we want to wait 60 minutes at this temperature. The following sequence will do it.



5. STEP "04." We want to ramp the temperature down from 90°C to 30°C in 90 minutes → ramp speed 45°/hour. The following buttons change the ramp speed to the new value.

ENT	8	4	5	ENT
-----	---	---	---	-----

6. STEP "05." We ramp down until 30°C with the following buttons.

ENT	ENT	1	3	0	ENT
-----	-----	---	---	---	-----

7. STEP "06." The bake profile is done, we stop the process with the special function HEATER OFF.

ENT	9	ENT
-----	---	-----

8. Buttons "7" and "8" can be used to check all the steps.
9. Exit the edition mode with button "6".