

Sawatec LSM-200 and HP-401-Z User Manual

Version of 2020-01-04.

1. Introduction

This user manual explains how to operate the Sawatec LSM-200 manual coater and the Sawatec HP-401-Z hotplate for the coating and drying/curing of SU-8 or polyimide (PI) films.

2. Login

- Login on "Z01 Sawatec LSM200" or "Z01 Sawatec HP401Z" with CAE on zone 01 accounting computer.


Z01 Sawatec HP401Z Hotplate for Soft-bake
Z01 Sawatec LSM200 Coater for negative resist

3. Sawatec LSM-200 coater operation

- First, users should edit the spin profile recipe according to the SU-8 or PI runcards.
- Recipe edition is accessed through the Sawatec touch-sensitive panel by pressing the following very intuitive ☺ sequence:



- Next, users will be able to edit the **Time [s]** and **Speed [rpm]** of all segments. Up/Down arrows are used to cycle through the segments.



- For the last step (usually segment N°8), make sure to activate **End Seg. YES** so that following segments are ignored.
- After edition, users can save (and later load) the recipe, with a name and description, by pressing this icon:




- After edition, users can go back in the menus with the "door" icon:



- Next, users will access the menu to start the spin-coating process with:




- Before going further, always check that the **"Arm" button is ON (greyed out)** so that the cover glass will cover the bowl during spin-coating.



- It is recommended to activate the "Vacuum" and check that the wafer chuck vacuum is close to -0.8 Bar.

- The SU-8 or PI resists are dispensed directly from the bottle. The resist pool should be approximately 5cm diameter.



- Make sure that the cover glass is not obstructed by any objects and start the process:



- Wait for process completion. The cover glass will move back to standby position. The wafer is then ready to be transferred to the hotplate.

4. Sawatec HP-401-Z hotplate operation


- Before transferring wafers to the hotplate, always double check that the actual temperature is not too high (from previous processed).
- Wafers are loaded on pins. Maximum 5 wafers are processed in the same batch.
- First, users should edit the softbake or PEB profile recipe according to the SU-8 (or PI) runcards.

Note: For PI baking, a “fast” softbake process at 65° and 105° is generally preferred. It is done on the free-to-use Süss hotplates close to the Sawatec units.

- Recipe edition is accessed through the Sawatec touch-sensitive panel by pressing the following very intuitive ☺ sequence:



- Next, users will be able to edit the **Time [s]** and **Temp. [°C]** of all segments. Up/Down arrows are used to cycle through the segments. The 1st segment will have both starting and final temperatures, and it will be generally used to do a RT relaxation before softbaking.



- For the last step (usually segment N°4 or N°5), make sure to activate **End Seg. YES** so that following segments are ignored.
- After edition, users can save (and later load) the recipe, with a name and description, by pressing this icon:



- After edition, users can go back in the menus with the “door” icon:




- Next, users will access the menu to start the spin-coating process with:



- Before going further, always check that the **“Heating” and “Cap” buttons are ON (greyed out)** so that the heating will be ON and the cap will close providing vapour confinement and exhaust during the bake.

In addition, make sure that the actual temperature is less than 5°C difference from the programmed starting temperature as the system will give an alarm.



- Start the bake process:



- Wait for process completion. The cap will open when done. The wafer is then ready to be stored in your box ready for the next step.

Note: The Sawatec hotplate does not have any active cooling so cooling down can take longer than programmed. It is safe to remove wafers when the temperature is < 70°C.

- Clean the working place properly.
- When done, logout with CAE on zone 01 accounting computer.