

Introduction to software

J.Ph.PIEL (PhD) Application Lab Manager SOPRALAB



II - WINSE SOFTWARE



II - WINELLI II SOFTWARE



III - SOPRA R&D SOFTWARE



EPFL – October 7th 2008



WINSE SOFTWARE

- Full measurements
- Basic analysis
- Recipe building

Two working modes :

- Engineer Mode
 - Buiding recipe
 - Measuement testing
 - Basic analysis testing
- Operator Mode
 - Loading recipe
 - Fully automatic measurement and analysis





Starting Screen





Main Screen





Main Screen



6



WinSE fo	or GES5E - ENGIR	IEER MODE - Product : default - Job : default *			
He Help	_				langed by a
	10/3/	2008 5:53:09 PM	DOPRA	011	Logged In : maintenance
MAST	TER			LOS	Maintenance SOPRA
		SYSTEM AWAITING PRO	CESS		Visible
Product :	Inknown	Job : default	E		Automatic Display Switch
	<u>N</u> ew recipe	ELLIPSOMETRY - Sample ref. Operation Analysis Sites Reports co	unfig		
	<u>O</u> pen	Coordinate Coird Dies C XY array		Descriptio	
	<u>S</u> ave	Water parameters Diameter (mm) C 100 C 150 C 200 C 300	Sample	Descriptio	
	Save <u>a</u> s	Fiducial type : C Flat C Notch			
	Save as default	Orientation angle (*) : () Sample thickness			
	New multiJo	Thickness (µm) : 700			
	Open				
	Save				
	Save as				
	Save as default				
Automatic	Recipe	Manual Result Setup Tools	Data Log Users Alarms		



WinSE for File Help	r GESSE - E	NGINEER MODE - Product : default - Job : default		
				ed In :
	1	0/3/2008 5:50:12 PM		enance
MAST	ER	SYSTEM AWAITING PROCESS	Log Maintenar	ICE SOPRA
			Vis	sible 4
Product : U	nknown	▼ Job: default ▼		Automatic Display Switch
	<u>N</u> ew reci	pe Sample ref. Operation. Analysis Sites Reports config.		
	<u>O</u> pen	MONOWAVELENGTH	peration	
	<u>S</u> ave	Incidence Angle Position Type Fixed Value (*)	Focus ○ Manual ○ Autofocus ○ None	asures library
	Save <u>a</u> s.	Fixed Fixed: 75	Autofocus per site	natif ox sca ms
	Save as <u>d</u> efault	○ Free	Site parameters Measures per site :	natif ox sca pb
	New mult	iJob Integration Time Type Fixed	Noise Measurement Never C Each Sample	ativeoxccdms
	Open	Fixed Variable Time(s) : 1 Position Type Scanned Values (µm) Unit	C Each Site C Each Acquisition C Periodic	nativeoxccdpb
	Save	CFixed Min 0.4 Max 0.8		hick ox ccd ms
	Save as.	C Scanned Step 0.1		hick ox ccd pb
	default	Fixed 45 Speed (Hz) : 6	Ontions	hick ox sca ms
		Use Mirror Analyser Tracking Mode	☐ Use MicroSpot ☐ Use Compensator	hick ox sca pb
		- None Current Previous	Correction Non Linearity Applied	hin ox ccd ms
			Save Measurement Automatically	thin ox ccd pb
Automatic	Recip	Manual Result Setup Tools Data Log Users	Alarms	



WinSE fo	or GESSE - ENGINEE	RMODE - Product : default - Job : defau				
	10/2/20	00 5-51-02 DM				Logged In :
	1013/20	56 5.51.65 P.M			tho post	maintenance Maintenance SOPRA
MINS	(LR	SYSTEM AV	VAITING PROCESS		100	Visible
Product :	Jnknown	▼ Job : <mark>default</mark>	Ē			Automatic Display Switch
B	<u>N</u> ew recipe	ELLIPSOMETRY Sample ref Operation Analysis	Beports config	Analy		
	<u>O</u> pen	Layers Regression				Regression Type
	<u>S</u> ave	Number of layers : 1		Definition	of the mo	Cel Cone Winse
	Save <u>a</u> s	Component Component	entration (%) Warning (%) Alarm (%) Min Max Min Max	Thickness (Å) Warning Value Min Mar	Alarm Min Max	Film Stacks library
	Save as <u>d</u> efault	1 Cauchy (n.k) PARAMETE 0		5000 0 0		reference-68
	New multiJob					
	Open					
	Save					
	Save as					
	Save as default	Substepto		Add new laver Delete	bottom laver	
		SILKIK		du new layer		
					Offsets	
Automatic	Recipe	Manual Result Setun	Tools Data Lon	lisers Alarms		

SOPRA





File Help	r ursor - triu	INTER MOUL - Product : default - Job ; defa				
	101	20000 5-50-24 DM	_			Logged In :
	100	3/2006 5.52.34 PM	-		tto pot	maintenance Maintenance SOPRA
Mints I		SYSTEM A	WAITING PROCESS		1	Visible
Product : U	Inknown	Job : default				Automatic Display Switch
T	New recipe	ELLIPSOMETRY		Sites de	scription	
	<u>O</u> pen	Sample ref. Operation Analysis Sample informations (mm)	Sites Reports config. Mapping informations			Mappings library
	Save	X origin 0.00 X incursion 57			X Y 1 10 10 2 10 380 2 10 380	25pts-10mm
	Save <u>a</u> s	Y origin 0.00 Y invension F Sites number: 25			4 10 1120 5 10 1490 6 467.5 10	
	Save as default	Mapping parameters	(++)	++++	7 467.5 380 8 467.5 750 9 467.5 1120	
	New multiJ	ob Vertical lines: 10 Horizontal lines : 8		+ + + + +	10 467.5 1490 11 925 10 12 925 380 13 925 750	
	Open			• + + + +	14 925 1120 15 925 1490 16 1382.5 10	
	Save				17 1382.5 380 18 1382.5 750	3
	Save as	C Radial		+ + + + + + + + + + + + + + + + + + +	19 1382.5 1120 20 1382.5 1490 21 1840 10 22 1840 20	
	Save as default			++++/	23 1840 380 23 1840 750 24 1840 1120 25 1840 1490	3
		Edge Exclusion (mm) Edge: 0 Flat: 0				
					Add Ins	ert
			Unselect All		Delete So	nt
Automatic	Recipe	Manual Result Setup	Tools Data Log	Users Alarms		

EPFL - October 7th 2008



Measurement Spectrum

File Help	R MODE - Product : default - Job : default *			
10/3/20	08 5:54:34 PM	JOPRA	011	Logged In : maintenance
MASTER	SYSTEM AWAITIN	G PROCESS	Log	Visible
Product : Unknown	▼ Job : default			Automatic Display Switch
Initialization	Status Focus Signal Measurement Graph Tan(Psi) - Cos(Delta) Intensity Power P	Results Graphics Power / Reference Alpha Beta Acquisition Data	Informati	on Configuration
		TAN (PSI)	Curren	nt Acquisition
			Incid	Wave 0.000
Start				ntensity 0.00 cts
	Tan(ps			Power 0.00 cts/s
			Anal	yzer Angle 0.00 *
		e.V		
		COS (DECIN)		
Abort				
	(etta)			
	Cos(c			
		۰V		
	MEASUREMENT IN PROGRESS			
Automatic Recipe	Manual Result Setup Too	Alarms		



WINELLI II SOFTWARE



For Simulation : Tan (Ψ), Cos (Δ), T_p, T_s, R_p,R_s For fitting using Levenberg Marquard approach

On Ellipsometric data: Tan (Ψ) and Cos (Δ) spectra On photometric data : T_p, T_s, R_p, R_s spectra

Using :

- Effective Medium approximation
- Polynomial law
- Harmonic oscillator function
- Drude Model
- Rough layer
- Anisotropic layer



WinElli_II ver 2.2.0.3



Starting Screen



EPFL - October 7th 2008



How to load a measurement

	🗲 Winelli II Software v2.2.0.8									
	Eile Ireeview Setyp Infos Help									
)			
	Session 1 SpreadSheet_1	_								
		<u> </u>	Sessi	sion 1 - Document			X			
		Eile	le E	📰 Session 1 - SpreadSheet_						×
Click with	right	(6	File Export Analysis GetValue		-				
buttor	٦			$\overline{\boldsymbol{\rho}}$		٢		(Å) (Å		
		6	765	Data list .		E	nergy (eV)	Tan Psi	Cos Delta	Alpha 🔺
			$\exists I$	Sample 1-2	-	1	0.7500	0.6211	0.4074	5.1260E-
Session 1			$/\Lambda$	Selected data tune :	Dimension -	2	0.7750	0.5962	0.4112	-4.0960E-
Add Function	Spreadsheet		A	Selected data type .		4	0.8250	0.5397	0.4393	-5.0150E-
Delete Session	Structure standard	Y /		Ellipsometry	2D	5	0.8500	0.5144	0.4403	-4.8040E-I
Delete Dession	N.O. K. Terdises Duilden			Commont :		6	0.8750	0.4883	0.4434	-5.1930E-
	N & K Indices Builder			comment.		7	0.9000	0.4689	0.4301	-4.0510E-
	Etching Calculation					9	0.9500	0.3941	0.3690	-9.9690E-
				Number of points - 17	T. Speetral unit : U	10	0.9750	0.3611	0.3360	-8.7340E-I
					T Specual unit. ev 👻	11	1.0000	0.3261	0.2827	-0.10
				Range : 0.750	0 to 5.000	12	1.025	0.2862	0.2353	-0.12
				,		14	1.030	0.2350	0.1730	-0.12
				Micro Spot : 🦵		15	1.100	0.1961	0.1773	-0.12
						16	1.125	0.1697	0.1811	-0.14
				- Measurement parameters for	r ellipsometru	17	1.150	0.1520	0.2429	-0.10
				Measurement parameters for		19	1.200	0.2075	0.3184	0.23
				Analyser Mode :	Current tracking	20	1.225	0.2593	0.3190	0.21
				Integration time Mode :	Fixed at : 5(s)	21	1.250	0.3123	0.3524	0.18
				Background Mode :	Meas he scan	22	1.275	0.3640	0.3839	0.15
				Nac Baratika ana Kina a	meds by sedn	23	1.300	0.4642	0.4187	0.10
	CIICK to open			Non inearly correction .	Hadamard correction	25	1.350	0.5140	0.4783	0.10
	· · ·			Others :	Nothing	26	1.375	0.5632	0.5014	9.1150E-
	one or several					27	1.400	0.6199	0.5294	9.5620E-1
				Parameters :		29	1.450	0.7360	0.5548	8.8870E-I
	data fila(a) in			Angle of incidence		30	1.475	0.8184	0.5630	0.10
				75.00 -		31	1.500	0.9220	0.5731	0.11
						32	1.525	1.067	0.5852	0.14
	current					34	1.575	1.601	0.6685	0.22
	o di i o i it					35	1.600	2.145	0.7970	0.28
	session					36	1.625	2.744	0.9591	0.24
	30331011			1						
						_	_			



How to display the measurement graph





How to display a default structure

Click on icons





How to define a model for each layer

🛫 Winelli II Software v2.2.0.8									
Eile Ireeview Setyp Infos Help									
	۲	6 🛞 🧐		(2))				
Session 1									
E SpreadSheet_1	Session 1 - S	tructure_2							
Structure_2	File Structure La	iyer Multilayers					-		
	-Structure param	neters for file : New Structure	•				-		
	Layer	Description	Thickness (um)	Optical index	Gradient Profile	Multilayers (*n)			
	Ambient	Void	G5	nk file 🗸	None 🗸				
	1	Dispersion la w	30.00	Dispersion la w 📿	None 🗸				
	Substrat	sio2		nk file Mixed nk files	None 👻				
				Dispersion law nk file & dispersion law					
				Alloy Diffusion layer					
				\wedge					
	Lawar att -								
	Optical index	type : Dispersion law	- Terms		- Find dis	enersion law			
	law tune · 🔽	tandard Dialaatria Eurotian	-	/					
	an grot ja	ranuaru preiectric Function	UVY	Modif Cauchy	/ [n.] 🗾 🚺	In Database			
	Nb of peaks :	0	19/	No model		Define law			
	5			\					
			7						
			/						
		/							
			Clic	ck to selec	t optica	I model :			
		Motorio	l of t	ha aurran	Elovor c	on ha de	ofined		
		_ materia	ιοιι	ne curren	l layer c	an be de	annea		
		by sing	le or	mixed N8	K files	tired fron	n l		
		by Sing							
		sopra's	data	a base or b	oy dispe	ersion lav	<i>w</i> s, allov		
			diff.		7		, ,		
		models	, anti	usion laye	IS				
her /ºº//UUX									



Dispersion Law builder





How to a add a Lorentz harmonic oscillator

Winelli II Software v2.2.0.8		
le <u>T</u> reeview Set <u>u</u> p Infos <u>H</u> elp		
	(2)	
Session 1 SpreadSheet_1 2D-Graphics_1 Structure_2	Session 1 - Structure_2 File Structure Layer Multilayers Dispersion Law Builder Ele Owner:	
	Dispersion law informations : Display : Display unit : Display unit : Display Range : Nb Points : Law type : Standard dielectric function	
	Comment :	
	Number of peaks : n* Model 1 Lorentz I Add Del I Display all Peaks I 1 Add 1 Display all Peaks I I 1 Add 0 1 2 Display all Peaks 1 A 1 A 2 Lambda0 0 0.5000 ymma I 3 Gamma 5.0000E-02, µm	



How to define a Drude Law

🗲 Winelli II Software v2.2.0.8		
File Ireeview Setup Infos Help		
	(2)	
Session 1 SpreadSheet_1 2D-Graphics_1 Structure_2	View View	
	0.2 0.4 0.6 0.8 1 1.2 1.4 1.6 1.8 Wavelength (m)	

How to enter in the regression parameters window





How to define the fitting parameters





Fitting procedure



EPFL – October 7th 2008



Result of the fitting procedure



EPFL – October 7th 2008

JOPRA

How to optimise measurements



Select the spectral range Select the spectral résolution Select the appropriate Angle Of Incidence





 \succ

How to optimise analysis

- Start with the simplest model
- Select the appropriate model for the optical index determination:

Bruggeman : Poly silicon, Roughness, Porous Materials

Dispersion law : need appropriate started values SDF-UV term and Lorentz peaks for dielectrics (most of the time) SDF-Drude and Lorentz peaks for metals Forouhi Model (FIM) for amorphous as SiONx Model Dielectric Function (MDF) for Alloy materials (SiGe, Si...)

Bulk calculation : for substrate or fully absorbing layers.

Point by point : need a well defined model (included T knowledgeEPFL - October 7th 2008and interfaces if necessary)



How to optimise analysis

Select the appropriate model for the regression
 Tan Psi&Cos Delta for « ultra thin layers »
 Alpha&Beta for « thin and thick layers »
 Pseudo Dielectric Functions: ε_r and ε_i useful for alloy material for instance.

Optimise the analysis

Adjust the spectral range

Add interface (rough layers) or/and gradient (inhomogeneous layers).

EPFL – October 7th 2008

Use special option (Aperture for microspot measurements...)



SOPRA R&D SOFTWARE



- Basic displacements of the different motor axis
- Basic acquisition of the signal : Show mode
- Basic measurements
- Calibration of the system
 - Spectrometer
 - Spectrograph
 - Non lineary of the detectors
 - Polariser and analyzer axis versus the plane of incidence : A0 and P0
- Hardware Setting



Starting Screen





Main Screen





Basic Displacements of each axis

SOPR	A R&D Software (Rev. 3.0.10) Administ	ator	
File Mea	surement Calibration Maintenance	Information	-	201 2010 9 14
	😵 Status: Initialized		🕤 Info	Info: End of measurement
D	Axis Tools		2	
9				
00				
	•		-	
			-	
	Init Send to	Cont Rotati	n	
	Polarise			
	Spectrome	er		
	Goniomet	er		
	Polariser A	m		
	Analyser A	m		
	Mapping		- 1	
	i) Spectrometer	r 68*	- 1	
	1) Spectrometer	r /6*	- 1	
	 Deterder 1 (Muel 	er Matrix)		
	 Retarder 2 (Muel 	er Matrix)	=	
	SOPRA			



Show Mode

File M	easurement Calibration Maintenance Information		
	Status: Initialized 📃 🕕	nfo: End of measurement	
		Sums Intensity Tan(Psi)/Cos(Delta) Alpha/Beta Options	
		10-Item# Intens	ity (cts/s) Symetry (%) Residual
→ 🕙	Sensor	9-	
	Mono C Multi C FT-IR C X-Ray	8-	
*	PM auto	7-	
	Grating auto Measurement type		
	Ellipsometry Photometry	ternsity	
	General options		
	MicroSpots Attenuator 1	3-	
	T Shutter 📕 Attenuator 2	2-	
	Retarder options © With & Without © With	1-	
	Correction	U – Hadamard Sums	
	Analyser tracking None Previous	10-	
	Correction	9- 9-	
	Special features		
	Four sums Integration time (s) O Eight sums Polariser Frequency		
	☐ Single shot ☐ Log data	■ enter 5-	
	✓ Background (cts.s-1)		
		3- 3-	
	Spectrometer	2- 2-	
	Set WaveLength	5 1- 1-	
			4 6 8
hic hu	utton	Acquisition number	Acquisition number
INS DU	IIION	Optical bench Mueller Matrix Mapping Tilt	
ent <mark>er i</mark> r	n the show mod	Move Analyser (*) 45 X Y Z Axis 4 Save position (*) 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0 0	lion
		Goniometer Move 7 Move 7 Move 7	
		Arms uncoupled Axis4 Axis4	2 Joystick
		Move XY XYZ XYZ Theta Position Posit	ing ion



Measurements Parameters Window



During Acquisition ...



EPFL – October 7th 2008



General Calibration menu

🖬 SOPRA R&D Software (Rev. 3.0.10) Administrator			
File Measurement Calibration	Maintenance Information	5.25	
Status: Initia	Calibration 🕨 🕨	Prism Grating Incidence Incidence X-Ray Attenuators X-Ray	
	Hardware SetUp SetUp Manager Stop & RE-Start Tool		
		Polariser OffSet Non linearity Ellipsometry Non linearity Photometry Attenuator Detarder (VIS_NID_ETID)	
		Retarder neutral axis (Mueller Matrix) Retarder parameters (Mueller Matrix) Polarisation extinction rate (FTIR)	
		Spectrograph response	



Polariser and Analyser offset calibration

🗔 SOPRA R&D S	oftware (Rev. 3.0.10) Administrator
File Measurement	Calibration Maintenance Information
💈 Sta	Reference sample A0/P0 Ctrl+A User Management Info:



Hardware setting

