

MEMO

Laboratory waste bins: Recommendations for solid contaminated waste

Version date	April 2022
Category	Waste bin
Keywords	Waste bin, solid contaminated waste, special waste

1 Purpose

The purpose of this communication is to summarize the rules and recommendations regarding waste bins for <u>solid contaminated waste</u>: gloves, paper, plastic lab consumables, etc. (excluding glass and sharp/pointy objects) contaminated with one or more of the following GHS pictograms.













This communication applies to <u>solid contaminated waste only</u>; it does <u>not apply to</u> biological or radioactive waste. Waste contaminated with *nanomaterials* must be double packaged (in a bag or a plastic container).

Waste contaminated with an oxidizer must be collected separately in an UN-certified container (use the white plastic buckets).

2 Rules & requirements

There is no universal waste bin for EPFL laboratories. Each research unit can purchase a model of its choice outside of Catalyse <u>AS LONG AS it fulfills the requirements listed below</u>. To be used with the official EPFL bags for solid contaminated waste only. See sections 3 and 4 for examples.

- 1. Waste bins made of metal are highly recommended, especially for chemistry laboratories.
- 2. Plastic bins are allowed.
- 3. Bag holders are allowed.
- **4.** Bins without lid <u>OR</u> a with a foot pedal opening lid.
- 5. Adapted to the size of the official EPFL waste bag for solid contaminated waste (90-110 L).



3 Examples and recommendations



Official EPFL waste bag for solid contaminated waste.

Note: these bags are gray in Lausanne and Neuchâtel, blue in Sion.



Metallic waste bin available at the Chemical stores on campus.



Metallic waste bin available through Catalyse (select Kaiser+Kraft catalog; ref# 521214).





Metallic waste bin without lid available through Catalyse (select Manutan catalog; ref# A497898).



Waste bag holder with a foot pedal opening lid available through Catalyse (select Manutan catalog; ref# A518379).

Ex: https://www.manutan.ch/fr/fab/poubelle-anti-feu.



4 Examples of waste bins NOT recommended or NOT allowed.



These yellow plastic bins should only be used for glass waste and or sharp/ pointy waste.

Note: No needle, no scalpel blade.



These carboard boxes are <u>ONLY allowed</u> for waste contaminated with genetically modified organisms. For BIO laboratories only.

Note: No glass, no needle, no scalpel blade.