

LBEM

Laboratory of Biological  
Electron Microscopy



# Cryo-Electron Microscopy method development for Life Sciences


Henning Stahlberg



May 11, 2021

# Short CV

## *Cryo-EM Method Development, Neurodegeneration*



3  
Henning Stahlberg

**Henning Stahlberg** (born 1965 in Germany. Swiss and German Nationality)

- 1993: Physics Diploma at the TU Berlin, Germany
- 1997: PhD in Cryo-EM of Membrane Proteins with **Dubochet** (UNIL) and **Vogel** (EPFL)
- 2002: Habilitation in Struct. Biology of Membrane Proteins with Andreas **Engel** (Biozentrum)
- 2003-2009: PATT & Assoc. Prof. of Molec. & Cell. Biology, UC Davis, CA, USA
- 2009-2021: Prof. of Struct. Biology, Biozentrum Basel and Director C-CINA, Uni Basel
- 2020– : Prof. of Physics at EPFL, Prof. Fund. Microbiol. at UNIL  
to set up the  
*Laboratory of Biological Electron Microscopy*  
and the *Dubochet Center for Imaging*

# Transmission Electron Microscopy: TEM

---

- Transmission imaging
- Electrons at 300 keV energy.
- 3D reconstruction through image processing.
- Resolution on proteins: 3 Å
- Sample: <300 nm thin






# A Resolution Revolution in cryo-EM from Direct Electron Detector Cameras




Virus  
frozen  
in ice

until 2012  
Photographed on a  
CCD camera or on film



Rotavirus, mockup images adapted from Grigorieff, *eLife*, 2013;2:e00573.

2013  
Photographed on a  
direct electron detector



since 2014  
*Filmed* on a  
direct electron detector






Image after  
alignment and averaging  
of 100 movie-frames.

Movement profile:




## Blobology



FMO membrane protein complex at ~2 nm resolution

## Atomic Structures



Apoferritin at 1.62Å resolution (Danev et al.)




# Cryo-EM Structure of a human ABC Transporter


ABCG2 with Fab fragments in lipid nanodisk

Data collection: 3 days on Titan in C-CINA,



ca. 3000 movies, ca. 200'000 particles, ca. 30'000 particles in final map




**a**




**b**




**c**



**d**





**e**



Taylor et al., Nature (2017), collaboration with Kaspar Locher (ETHZ)

# LBEM & DCI



1/2022

5/2021



**Talos (200kV)**  
With STEM



**Prototype Titan Krios (300kV)**  
E-CFEG, RF-Cavity, Cs-probe corrector, Falcon4

5/2021



**Titan Krios (300kV)**  
E-CFEG, SelectrisX filter, Falcon4

6/2021



**Titan Krios (300kV)**  
E-CFEG, Falcon4

6/2021



**Glacios (200kV)**  
X-FEG, Falcon4




# The LBEM Instrument: RF Cavity: Producing single electrons at GHz rate

Titan Krios G4 300kV, CFEG, F4 camera, RF Cavity, Probe Cs corrector



*Single electron illumination reduces the beam damage 2.5x or more (?)*




# The LBEM Instrument: RF Cavity: Producing single electrons at GHz rate

Titan Krios G4 300kV, CFEG, F4 camera,  
RF Cavity,  
Probe Cs corrector



*Single electron illumination reduces the beam damage 2.5x or more (?)*



# The LBEM Instrument: RF Cavity: Producing single electrons at GHz rate

Titan Krios G4 300kV, CFEG, F4 camera, RF Cavity, Probe Cs corrector



Single electron illumination reduces the beam damage 2.5x or more (?)

## Paraffin

3 GHz frequency





Figure 3. Mitigation of radiation damage at different duty factors of pulsed electrons. (a) Selected continuous- and pulsed-beam diffraction patterns of a paraffin crystal with the 110 family Bragg peaks highlighted. (b) Fading curves from the 110 family reflections and exponential decay fits as a function of pulsed beam duty cycle factor at 5.2 GHz repetition rate for the same dose rate (0.2546 e<sup>-</sup>/Å<sup>2</sup>/frame), grey guideline indicating critical dose ( $D_{crit}$ ).

## Bacteriorhodopsin




Figure 4. Comparison of the continuous random and pulsed random electron beam damage for the same dose rate (0.1045 e<sup>-</sup>/Å<sup>2</sup>/frame). (a) Selected continuous- and pulsed-beam powder diffraction patterns of a purple membrane with the (430) rings at total accumulated dose 0.014 e<sup>-</sup>/nm<sup>2</sup> (1.4 e<sup>-</sup>/Å<sup>2</sup>). Rings are diffuse, possibly due to sample thickness or vibrational blurring from soft phonon modes at room temperature. (b) Fading curves from the (430) powder rings and non-linear exponential decay fits for each measured duty cycle factor from 10 to 100%.

Strobo 300 Pulser  
Euclid Techlabs,  
Bolingbrook, Illinois, USA

Lau et al., (2020)  
Rev. Sci. Instrum. 91, 021301





# The LBEM Instrument: Coherent Electron Diffractive Imaging

Titan Krios G4  
300kV, CFEG, F4 camera,  
RF Cavity,  
Probe Cs corrector




*Electron diffraction  
with one protein at a time*

## Conventional electron source



# The LBEM Instrument: Coherent Electron Diffractive Imaging

Titan Krios G4 300kV, CFEG, F4 camera,  
RF Cavity,  
Probe Cs corrector




With  
Radosav Pantelic, EPFL,  
Bernd Schmitt, PSI


**Goal: Diffractive Imaging to obtain *phase-contrasted* diffraction patterns, which give more SNR and no resolution loss from drift.**

# The LBEM Instrument: Coherent Electron Diffractive Imaging

Titan Krios G4 300kV, CFEG, F4 camera,  
RF Cavity,  
Probe Cs corrector





The sample is exposed at multiple locations with the electron beam:



With  
Radosav Pantelic, EPFL,  
Bernd Schmitt, PSI



# LBEM & DCI




1/2022


5/2021



**Talos (200kV)**  
With STEM



**Prototype Titan Krios (300kV)**  
E-CFEG, RF-Cavity, Cs-probe corrector, Falcon4



5/2021



**Titan Krios (300kV)**  
E-CFEG, SelectrisX filter, Falcon4

6/2021



**Titan Krios (300kV)**  
E-CFEG, Falcon4

6/2021



**Glacios (200kV)**  
X-FEG, Falcon4





## Cubotron Ground Floor

Big thanks to Olivier Schneider (IPHYS) for making space available!



# Acknowledgements

Pratibha Kumari,  
Roland Riek, ETHZ

Paola Picotti, ETHZ

Beat Meier, ETHZ  
Anja Böckmann, Univ. Lyon

Ronald Melki, Luc Bousset, CNRS Paris

Lawrie Rajendran, UniZH

Marialuisa Quadri,  
Wilfred van Ijcken,  
Tim Moors,  
Evelien Huisman,  
Angela Ingrassia  
Yvonne de Gier,  
Vincenzo Bonifati,  
Annemieke J.M. Rozemuller,  
Wilma van de Berg, VUMC, Amsterdam

Pierre de Rossi  
Magdalini Polymenidou, Uni ZH

Anne De Paepe,  
Johannes Erny,  
Daniel Mona,  
Andreas Staempfli,  
Jörg Hörschemeyer,  
Gregor Dernick,  
Bernd Bohrmann,  
Matthias Lauer,  
Markus Britschgi, Hoffmann-La Roche, Basel

Nicholas Taylor, Novo Nordisk, Copenhagen  
Ricardo Guerrero-Ferreira, Emory Univ., Atlanta  
Sarah Shahmoradian, PSI, Villigen

Alexandra Graff,  
Christel Genoud, FMI Basel

Frederik Grossrüschkamp,  
Daniel Niedieker,  
Samir F. El-Mashtoly,  
Klaus Gerwert, Ruhr University, Bochum

Sabine Ipsen,  
Jürgen Hench,  
Gabriel Schweighauser,  
Stefan Frank, Pathology, Uni Basel  
Markus Tolnay, Pathology, Uni Basel

Julia Kowal  
Kaspar Locher, ETHZ

Karl-Heinz Altmann, ETHZ

Erik Fröjd  
Bernd Schmitt, SLS Detector Group, PSI

Jan Pieter Abrahams, Uni Basel and PSI

Andrea Ablasser, SV  
Pierre Gönczy, SV  
Hilal Lashuel, SV  
Suliana Manley, IPHYS  
Fabrizio Carbone, IPHYS  
Michael Unser, STI

**Funding: SNI, NCCR TransCure, SNF,  
Synopsis Foundation, Heidi Seiler Foundation,  
Hoffmann La-Roche, SystemsX.ch, Loterie Romande**

Carolin Böing  
Domenic Burger  
Marta di Fabrizio  
Amanda Lewis  
Inay Mohamed  
Dongchun Ni  
Radosav Pantelic  
Jason Selby  
Daniel Stähli

Mireille Fasmeyer

Sandor Kasas  
Sergey Sekatski  
Alexander Verkhovsky  
Kamilla Guérin  
Anne-Celine Kohler  
Zeno Messi  
Maria Inés Villalba  
Anton Malovichko



Heidi Seiler  
Foundation

Stiftung Synapsis  
Alzheimer Forschung  
Schweiz

